medinova

F12

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EN WARNINGS

HEALTH RISKS

Can cause irritation during contact with skin and eyes. Prolonged exposure to the product may cause allergic reactions.

PRECAUTIONS FOR USE

- Do not use if the product or the packaging is damaged, or if the appearance of the product has changed.
- Do not use after the expiry date.
- Respect the dilution ratio of the bag of F12, to avoid altering the effectiveness of the extender.
- Do not store near food, beverages or pet food.
- Wear gloves and a protective mask when handling the powdered product.

STORAGE

The product can be transported and stored at room temperature not exceeding +20°C. Do not freeze the powder or the diluted product.

PREPARATION INSTRUCTIONS

- To obtain the final solution, it is necessary to dilute the entire content of the bag of F12 in sterile double-distilled water type II according to the size of the bag (1-5-60-100 lt).

- It is important to accurately measure the volume of sterile bidistilled type II water to avoid any variation in osmotic pressure.

USE

- Dilute the entire content of the package in double-distilled type II water at 37°C.

- After verifying the completed dissolution of the powder, allow the solution equilibration for at least 30 minutes before adding the semen.

- Check the temperature of the sperm and solution to avoid thermal shock.

- Do not mix the medium and semen if the temperature difference between the two is greater than 2°C.

- Use the diluted medium in a maximum of 48 hours after preparation.

- After preparation of the medium, it is recommended that the unused solution be refrigerated between +2° and +6°C in a sterile, sealed recipient.

- Make the dilution of the semen by adding the necessary medium to obtain the desired concentration of live spermatozoa per 100 ml volume, depending on the quality and volume of the ejaculate.

- Store the doses between (+15°C / +17°C) in a thermostatically controlled refrigerator.

- After dilution use the doses within 12 days.

QUALITY CONTROLS

Physico-chemical characteristics of the product diluted in sterile bidistilled water type II.

- Osmolarity: 290 +/- 20 mOsm/kg H₂O

- pH: 7.2 +/- 0.2

PRODUCTION STANDARDS

The semen extender is manufactured in MEDI NOVA laboratories, an ISO 9001 2015 certified company. The extender complies with all quality standards, respecting the specifications established by European legislation.

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