# medinova

# F12

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## EN WARNINGS

### INDICATIONS FOR PREPARATION AND USE

- Dilute the entire contents of the bag of F12 in sterile double-distilled type II water at 37°C depending on the size of the bag (1-5-60-100 lt).
- Measure accurately the volume of type II sterile double-distilled water to avoid any variation in osmotic pressure.
- Check the complete dissolution of the powder.
- Let the solution equilibrate for at least 1 hour in a water bath at 37°C before adding the sperm.
- Do not mix the diluent and the sperm if the temperature difference between the two is greater than 2°C.
- Calculate the optimal concentration of spermatozoa/ml in the dose and add the diluent to obtain the final volume.
- Store the dose at (+15°C / +17°C) in a thermostatic fridge.
- After dilution, use the dose within 12 days.

### PHYSICAL-CHEMICAL CHARACTERISTICS OF THE PRODUCT

The physical-chemical characteristics of the product diluted in type II sterile double-distilled water at 37°C, after 1 hour are as follows:

- Osmolarity: 290 +/- 20 mOsm/kg H<sub>2</sub>O
- pH: 7.02 +/- 0.2

### PRECAUTIONS FOR USE

- Do not use in case of deterioration of the product, of its packaging or if the appearance of the powder has changed.
- Do not use after the expiry date.
- Wear gloves and a protective mask when handling the powdered product.

### CONSERVATION

- The product can be transported and stored at room temperature.
- Do not freeze the powder or the diluted product.
- After preparing the diluent, it is recommended that the unused solution be refrigerated between +2°C and +6°C in a sterile, sealed container and used within 48 hours.

### HEALTH RISK

The product can cause irritation in contact with the skin and eyes. Prolonged exposure to the product can cause allergic reactions. Do not store near food, drink or animal feed.

### PRODUCTION STANDARD

This product is manufactured in the laboratory of MEDI NOVA, a company certified ISO 9001 2015. The product conforms to all quality standards, respecting the specifications established by European regulations.

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